

Introduction to metagenome-based pathogen identification

Sahar Abdelrazek, M.Sc., Ph.D.


Department of Biomedical Sciences and Pathobiology

VA-MD College of Veterinary Medicine


Virginia Tech



Genome/metagenome-based pathogen identification




Methods for Pathogen Identification




Culture-based Methods

Time-consuming
some bacteria are difficult
to culture



Molecular Methods

Only pathogens with
known sequence data can
be detected




Sequencing Methods

Can detect newly
emerging pathogens as
well as non-culturable
pathogens

Metagenome-based pathogen identification


Metagenomics is the sequencing of a whole community of microorganism DNA as opposed to sequencing just the individual microbes.



Reasons for choosing metagenomics

- No need of culturing the pathogen
- Higher resolution
- Can detect pathogens from newly emerging diseases

Metagenome-based pathogen identification pipeline



1 Sample Collection

Planning & Preparation

Define Objectives: Study purpose (e.g., disease diagnosis, pathogen diversity, epidemiology).

Select Plant and Tissue Type: Choose specific tissues (e.g., leaves, roots, stems) based on where the pathogen is likely to be present.

Sample Collection

Timing: Collect samples at the right stage of disease development.

Tools and Materials: Sterile gloves, scissors or scalpels and sterile bags or tubes.

Storage & Transport

Transport: Ensure samples remain cool during transport to the lab (ice packs or a portable refrigerator if possible).

Immediate Storage: (-80)



2 | DNA extraction



community.nanoporetech.com/docs/prepare/extraction_protocols

- 25 sample types
- 72 protocols
- DNA and RNA
- Sample handling & storage advice

Role # 1: You can't sequence what you don't have?

Sample Preparation: In the lab, further clean the samples if needed. Use sterile conditions to prevent contamination.

Extraction Kit: Use a DNA extraction kit suitable for your sample. (modifications to handle each type of sample).


Purity Check: Use Nanodrop to check the DNA purity (A260/A280 ratio should be ~1.8).

Quantification: Use Qubit for accurate DNA quantification.

Sample
Preparation


Extraction Kit

Quality & Quantity
Assessment



3 | Library preparation

Role # 2: Trade-off between sequencing output and read length



Long Reads
Portability
Real-Time Data
Versatility



Any fragment length
libraries 20 bases to 4
Mb+




Rapid prep in just 10
minutes



PCR-free protocols to
eliminate bias

	Output optimised	Speed optimised	Ultra-long reads optimised
	Ligation Sequencing Kit	Rapid Sequencing Kit	Ultra-long DNA Sequencing Kit
Preparation time	60 minutes	10 minutes	200+ mins + 1x O/N incubation
Input recommendation	~1000 ng gDNA or 100-200 fmol for amplicons	~100 ng gDNA	6M cells
Fragmentation	Optional	Transposase-based	Transposase-based
Amplification	No	No	No
Barcode options	Native Barcoding Kit 24 Native Barcoding Kit 96	Rapid Barcoding Kit 24 Rapid Barcoding Kit 96	-
Typical output	● ● ●	● ● ○	● ● ○
Adaptive sampling	Yes	Yes	Yes


Protocol builder



Provides a range of end to end workflows to help select the best library prep solution


Which DNA library prep kit do I choose?

Let our tools help you choose the best solution




The Protocol Builder provides a range of end-to-end workflows to help select the best library prep solution:
community.nanoporetech.com/knowledge/protocol_builder


3 Library preparation workflow



Native Barcoding Kit 24 V14



3 | Sequencing



3 Real-time data analysis



EPI2ME

EPI2ME WIMP rapid species
identification and quantification
from metagenomic samples



3 | Example for real-time data analysis using WIMP

Tomato sample with wilt symptoms





1 | Sample
Collection

2 | DNA
extraction

3 | Sequencing &
Data analysis


4 | Pathogen
identification

Sample ID	22-263
Collect Date	October-12-2022
Virginia District	Southwestern
County	Buckingham, Virginia
Host Common Name	Tomato
Variety	Red Morning
Host Genus	Solanum
Host Species	lycopersicum
Oxford Nanopore Technologies	PromethION
Samples run per flow cell	20



3

WIMP Results



READS ANALYSED
1,088,000

READS CLASSIFIED
520,840

READS UNCLASSIFIED
567,160



Tomato sample with wilt symptoms

SHOWING TAXA WITH READS
NCBI Taxonomy Tree

MINIMUM ABUNDANCE CUTOFF


3% 1% 0.5% 0.1% 0%

4 Pathogen identification

HIDE FILTERS

SHOW TOP N TAXA

10 20 30 100 200



Taxon	Cumulative Reads
Ralstonia solanacearum	288,732
Homo sapiens	48,735
Stenotrophomonas maltophilia	30,867
Bordetella petrii	20,110
Enterobacter sp. ODB01	14,977
Acinetobacter baumannii	11,183
Agrobacterium tumefaciens	9,931
Enterobacter kobei	8,891
Pseudomonas aeruginosa	4,436
Clostridium tetani	3,227
Enterobacter cloacae	3,089

3

Run Kraken2 using PlusPF Kraken2 database

Classified species within our sample were then sorted by z-score, with the species having the highest z-score considered as the pathogen

Extract Ralstonia reads using Kraken

Assemble pathogen genome using Flye



genomeRiv

Read-based identification still

Genome-based identification

46 Archaea

Sankey diagrams showing the absolute microbiome read counts within our infected sample



Circular map of *R. solanacearum* assembled draft genomes visualized using GenoVi.



Genome/metagenome-based pathogen identification



GenomeRxiv for deeper genome taxonomic identification




Branchwater for large-scale sequence search and profiling



Core genome tree Phylogenetic analysis
AMRFinderPlus for Virulence and AMR prediction




LIntax for read-based strain-level pathogen detection



Explore



(Meta)Genomics-Based Pathogen Identification Center for Animal and Plant Disease Diagnostics, Biosecurity, and Pandemic Prevention



Our Team




**Kevin Lahmers, DVM,
PhD, DACVP—Anatomic
Pathology**

Clinical Professor of
Anatomic Pathology
205 Duck Pond Drive
Blacksburg, VA 24061
klahmers@vt.edu



Boris Vinatzer, PhD
Professor of Plant Pathology

551 Latham Hall
Blacksburg, VA 24061
vinatzer@vt.edu



**Lina Rodriguez
Salamanca, PhD**

Instructor and Plant Disease
Clinic Manager and
Diagnostician
114 Price Hall
170 Drillfield Drive
Blacksburg, VA 24061
linar@vt.edu



Sahar Abdelrazek, PhD

Research Assistant Professor
and
Project Manager
205 Duck Pond Drive
Blacksburg, VA 24061
abdelrazek@vt.edu

Available Services





**Comprehensive and
Rapid Pathogen
Screening**

**Custom Metagenomic
Analysis**



Research & Collaboration



Next: **genomeRxiv**

Questions

